

Shoulder Replacement Systems



Contact Info

- +971-50-9401025
- info@iltmedical.com
- www.iltmedical.com



Scan the QR code for
our full product catalog.




Reverse Total Shoulder Solution

V A I O S

SHOULDER SYSTEM
7A ODEP Rating

**Minimum
inventory**

Minimum components reduced complexity
without compromising clinical choice.




Cortical Bone Screw
±15° Adjustable
Superior/Inferior
Locking Screws



Glenoid Dome
Polished Cobalt
Chrome (CoCr)



Inverse Humeral Insert
Ultra High Molecular
Weight Polyethylene
(UHMWPE)



Locking Screw

Humeral Metal Back
Titanium



Humeral Neck
H-A.C. Zoned Coated Titanium
Fixed angle – 140° Neck Stem Angle
(NSA)



Humeral Stem H-A.C.
H-A.C. Supravit® Coated Titanium
8, 10 and 12, 14 and 16mm Diameters

